Modelo BBM-100

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Model BBM-100

Instruction

Thank your purchase of BBM-100, the Hand Massager. Before use, please read the manual instruction thoroughly, especially the caution part, so that you will operate it correctly, which should be well kept for your reference at any time.

Contents

The features and functions with usage

Health and cautions

Suggested usage

Maintenance and storage

Specification

The features and functions with usage

- 1 .Please insert the plug into the electric socket;
2. Push the control switch to adjust the speed and intensity; up for fast, down for slow;
3. Adjust the speed and intensity against your interest, enjoy the massaging;
4. The system of average weight will flap and massage equality to achieve saving time and energy;
5. The apparatus with light weight and strong functions can massage the parts where we are difficult to be;
- 6.The special designed mat of beating, offer the central massaging effect, make the massaging more comfortable and relaxing;
7. It is quiet when operating.
8. Infrared acupuncture will effective penetrate to correspond point, and strengthen the stimulation.
9. Please operate the infra red function by the switch in the underside of unit.
10. The usage of replacing massage- pinpoint. Press Power to shut (OFF), unplug from the socket, then screw the intrinsic massage-pinpoint off by counter-clockwise entirely, and screw up the needed pinpoint on the massager by clockwise, now you can use it.

Health and cautions

- ~ After use and before clean it, please turn it off and pull out the plug;
- ~ Please do not sit or stand on it hard;
- ~ Don't use when make shower;
- ~ Don't immerse any part of the massaging in eater;
- ~ Do not us near to an oven or any other heating element;
- ~ If the massager was dipped into water, you should pull out the plug immediately, do not take it out from the water directly;
- ~ Please operate it in the condition of circulating air, do not use it under the coverlet, pillow ,otherwise the machine may become too hot, causing firing, hidden problems of fire or hurt the user;
- ~ Do not use it on the hypersensitive skin;
- ~ The round do not touch your head, chest knee and elbow joint directly, to avoid body's hurt;
- ~ Use the massager only after resting for one hour after each meal, and drink 150cc water prior to massaging and 250cc after it;
- ~ It is advisable to use the massager for only 30 minutes each time;
- ~ Hold the socket when unplugging, do not pull the wire;

Anyone undergoing medical treatment or suffering from the following conditions should consult a physician before using the massager as instructed:

1. Expecting mothers ,customary abortions and unhealed or inflamed post-surgical wounds;
2. Acute diseases, heart diseases, bone fractures, dislocated joints, or individuals suffering from physical anomalies;
3. The feet with turgescence and bleb

Cautions:

Please don't use again until it is stopped 20 minutes, to avoid the machine becomes too hot causing the matters.

Suggested Usage

Neck/Shoulder

Thighs/Calives

Waist/Back

Feet

Maintenance and storage

- ~ Make sure that the plug had removed before cleaning;
- ~ Never attempt to use water, thinner, benzene or alcohol to clean the unit;
- ~ For cleanliness, use a damp cloth or towel to clean it, the wet one is prohibited;
- ~ The massager equipment is designed specially, no need to lubricate and maintenance in order;
- ~ If you will not use it in some days, please pack the wire and unit on the place of dry and far from fire
- ~ Please place it far from the heat and fire, do not place it under the sun long time.

Specifications:

Description: Multifunctional Hand Massager

Type: BBM-100

Voltage: 220V/50HZ

Frequency: 3200rpm (fast) 2800rpm (slow)

Power: 0.15A/28W